trovaprezzi.it

beurer FT 09/1

DEUTSCH

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ENGLISH

Signs and symbols
The following symbols are used in these instructions for use on the packaging and on the type plate for the device and accessories:

Caution
Follow instructions for use

Type BF applied part
Disposal in accordance with EC Directive WEEE (Waste Electrical and Electronic Equipment)

Manufacturer
Dispose of packaging in an environmentally friendly manner

Protected against foreign objects equal to or greater than 12.5 mm and against temporary immersion in water

Batch designation
Permissible storage and transport temperature and humidity

Operating
Permissible operating temperature and humidity

The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical devices.

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– part 3: Performance of compact electrical thermometers (non-predicative and predictive) with maximum device and the European Standard EN60601-1-2 (in accordance with CDRP 11, IEC 61000-4-2, IEC 61000-4-3 and IEC 61000-4-8), and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this unit.

Important safety information – retain for later use.
In order to make optimum use of all the merits of the clinical thermometer, you should carefully read the operating instructions before using it, retain the instructions for later use and also keep them accessible for other users.

The thermometer is intended to measure the temperature of the human body.
The thermometer is only designed for the measuring area on the human body stated in the instructions for use.

Children should not use the thermometer without supervision.
Before each use, check the thermometer for signs of damage or wear. Damaged or worn thermometers should not be used.

Always wait for the acoustic signal which indicates the minimum time required for measurement. Discuss the temperature with your family physician.

The thermometer contains sensitive electronic parts. Therefore, protect it from blows, bending, high temperature or direct exposure to sunlight.

Using the thermometer in strong electromagnetic fields such as next to a mobile phone, can cause a malfunction.
When the thermometer is switched on, it carries out a self-test. Therefore the accuracy of measurement need not be tested at regular intervals.

If the temperature is very abnormal, please contact your family doctor promptly.
This device meets the requirements of the EU Directive 93/42/EEC concerning medical devices, as well as those of the Medizinproduktegesetz (German Medical Devices Act), the European Standard EN 12470-3: Clinical Thermometers

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Measurement in the mouth cavity (oral)
For this purpose, guide the tip of the thermometer carefully into one of the heat pouches beneath the tongue, to the left or the right of the root of the tongue.

Measurement of temperature in the armpit (axillary)
This method of measurement is relatively inaccurate. Therefore, it cannot be recommended from a medical point of view.

Mode of measurement
Recommended time required for the measurement
Mean deviation to the rectal method

in the anus until the acoustic signal
in the oral cavity until the acoustic signal
in the armpit minimum 5 minutes

Use
First remove the protective foil from the display. For switching on, briefly press the on/off key. A short acoustic signal will announce that the thermometer is switched on. The thermometer will carry out a self-test for about 2 seconds. During the self-test all segments of the display will appear. Then the reference reading 37°C appears. Then the measurement symbol "C" and "Lo" is shown on the display. The thermometer is now ready for use. During the measurement, the current temperature is displayed continuously, and the "C" symbol flashes on and off. The measurement is completed when a stable temperature has been reached. When this signal has been reached, an acoustic signal is heard, the "C" symbol stops flashing and the temperature measure is displayed. As long as the measured temperature is below 37°C, "Lo" will appear. Higher temperatures will be shown directly, e.g. 32.1°C. By prolonging the length of measurement until after the acoustic signal, more accurate measurements can be obtained. In case the temperature is measured in the armpit, the thermometer should be left in the recommended

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