Booklet, 148 x 210 mm, 30 pages (6 pages cover), 2C, black + Pantone 300u

P&G RESTRICTED

# PROFESSIONAL CARE

designed by **BRAUN**

## OXYJET

Deutsch	4
English	7
Français	9
Español	12
Português	14
Italiano	16
Nederlands	19
Dansk	22
Norsk	24
Svenska	26
Suomi	28

<b>DE</b>	<b>AT</b>	<b>00 800 27 28 64 63</b>
<b>CH</b>		<b>08 44 - 88 40 10</b>
<b>GB</b>		<b>0 800 731 1792</b>
<b>IE</b>		<b>1 800 509 448</b>
<b>FR</b>		<b>0 800 944 803</b> (service et appel gratuits)
<b>BE</b>		<b>0 800 14 592</b>
<b>LU</b>		<b>800 21172</b>
<b>ES</b>		<b>900 814 208</b>
<b>PT</b>		<b>808 20 00 33</b>
<b>IT</b>		<b>800 124 600</b>
<b>NL</b>		<b>0800-4 45 53 88</b>
<b>DK</b>		<b>70 15 00 13</b>
<b>NO</b>		<b>22 63 00 93</b>
<b>SE</b>		<b>020 - 21 33 21</b>
<b>FI</b>		<b>020 377 877</b>

**Internet:**

[www.oralb.com](http://www.oralb.com)  
[www.braun.com](http://www.braun.com)  
[www.service.oralb.com](http://www.service.oralb.com)

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# English

## Welcome to Oral-B!

Before operating your Oral irrigator, please read these instructions and save this manual for future reference.

## IMPORTANT

- The installation instructions shall state that parts that have to be fixed must be fixed so that they cannot fall into water, unless they are of IPX7 construction.
- Periodically check the entire product/cord for damage. A damaged or non-functioning unit should no longer be used. Do not modify or repair the product. This may cause fire, electric shock or injury.
- If the product/ cord is damaged, take it to an Oral-B Service Centre.
- Oral irrigators can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge under supervision.
- Oral irrigators shall be used only with liquids defined in the manufacturer's instructions.
- Do not use iodine, saline (salt solution), or water insoluble concentrated essential oils in the appliance. Use of these can reduce product performance and will shorten the life of the product. Use of these products may void warranty.
- Use this product only for its intended use as described in this manual. Do not use attachments which are not recommended by the manufacturer.

## WARNING

- When unplugging, always pull the power plug instead of the cord. Do not touch the power plug with wet hands. This can cause an electric shock.
- If you are undergoing treatment for any oral care condition, consult your dental professional prior to use.
- This irrigator is a personal care device and is not intended for use on multiple patients in a dental practice or institution.
- Do not chew or bite on nozzle. Do not use as a toy or without adult supervision.

- For water flossing only.
- Do not direct water into the nose or ear. The potentially deadly amoeba, *Naegleria fowleri*, may be present in some tap water or unchlorinated well water and may be fatal if directed into these areas.

## Description

- a Nozzle
- b Nozzle switch
- c Nozzle compartment
- d Nozzle release button
- e On/off switch
- f Water intensity control dial
- g Handle switch (water stop and go)
- h Irrigator handle
- i Water container (600 ml)
- j Cord storage
- k Wall holder

## Specifications

For voltage specifications refer to the bottom of the charging unit.

Noise level: 76 db (A)

## Using your irrigator

Your Oral irrigator comes with a convenient wall holder. It can be mounted easily using the attached items (picture 1).

Before using for the first time, fill the water container (i) and allow the water to run through completely by holding the irrigator handle (h) over a washbasin.

1. Attach one of the irrigator nozzles (a) to the handle until it locks into place (picture 2). Fill the water container with luke warm water.
2. Turn on the irrigator by pressing the on/off switch (e) (picture 3). (When doing so, the handle switch should be in the lower position = water stop).
3. Select the water pressure with the pressure control dial (f) (1 = gentle, 5 = strong) (picture 4). We recommend that you start with setting 1, especially if your gums tend to bleed or if children are using the Oral irrigator.
4. Bend over a washbasin and guide the nozzle to your teeth and gums. Slide the handle switch (g) up to start the water stream (picture 5). Keep your mouth slightly open to allow the water to run into the washbasin.

Slide the handle switch up (g) for a continuous water stream.

Slide the handle switch down and press into it, to use the on-demand function (picture 6).

**Caution:** Turn off the Oral irrigator by switching off the on/off switch once the water container is empty. Longer operation of the Oral irrigator without water may damage the appliance.

## OxyJet nozzle

### Nozzle functions:

You can choose between two different nozzle functions by sliding the nozzle switch (b) down or up. Before doing so, turn off the water stream by sliding the handle switch down:

Rotating stream for general cleaning and gum massage (nozzle switch in lower position): The rotation of the water stream and its micro air-bubbles are designed for easy penetration into the interdental space and gentle gum massage (picture 8).

Focused stream type (nozzle switch in upper position): The focused water stream with its micro-bubbles are designed for targeted cleaning (picture 9).

## WaterJet nozzle

### Jet functions:

You can choose between different stream types by click-in and click-out the blue center of the nozzle. Before doing so, turn off the water stream by sliding the handle switch down:

When the blue center is out (default position) you will get a high pressure single precision stream (picture 10).

As soon as you click the blue center in, you will get a multi-jet stream, designed for great reach and massage effect of the gums and the gumline (picture 11).

Maximal operation time: 15 minutes

Cooling off time: 2 hours

## After use

Always empty the water container entirely. Keep the Oral irrigator switched on until the water stops running through the nozzle in order to prevent accumulation of bacteria. As soon as no more water flows out of the nozzle, switch the device off (picture 3). Press the nozzle release button (d) (picture 7) to take off the nozzle from the Oral irrigator. Always wipe the handle and the Oral Irrigator dry; place the handle back onto its base and store the Oral irrigator in the nozzle compartment (c).

## Cleaning recommendation

After use, rinse the nozzle under running water. Regularly clean the housing with a damp cloth. The water container and the nozzle compartment are removable and dishwasher-safe.

Subject to change without notice.

## Environmental Notice

Product contains batteries and/or recyclable electric waste. For environment protection do not dispose in household waste, but for recycling take to electric waste collection points provided in your country.

### Guarantee

We grant 2 years guarantee on the product commencing on the date of purchase.

Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose. Depending on availability replacement of a unit might result in a different colour or an equivalent model.

This guarantee extends to every country where this appliance is supplied by Braun / Oral-B or its appointed distributor. This guarantee does not cover: damage due to improper use, normal wear or use, especially regarding the nozzles, as well as defects that have a negligible effect on the value or operation of the appliance.

The guarantee becomes void if repairs are undertaken by unauthorized persons and if original Braun / Oral-B parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorized Oral-B Braun Customer Service Centre.

This guarantee in no way affects your rights under legal/ statutory law.

### Warranty replacement nozzles

Oral-B warranty will be voided if irrigator damage is found to be attributed to the use of non-Oral-B replacement nozzles.

Oral-B does not recommend the use of non-Oral-B replacement nozzles.

- Oral-B has no control over the quality of non-Oral-B replacement nozzles. Therefore, we cannot ensure the cleaning performance of non-Oral-B replacement nozzles, as communicated with the irrigator at time of initial purchase.
- Oral-B cannot ensure a good fit of non-Oral-B replacement nozzles.
- Oral-B cannot predict the long-term effect of non-Oral-B replacement nozzles on irrigator wear.

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