L I F E I N B A L A N C E

SOEHNLE

BEDIENUNGSANLEITUNG	2
USER MANUAL	4
MODE D'EMPLOI	6
ISTRUZIONI PER L'USO	8
GEbruIKSAANWIJZING	10
INSTRUCCIONES DE MANEJO	12
MANUAL DE INSTRUÇÕES	14
BRUKSANVISNING	16
BRUGSANVISNING	18
INSTRUKCJA OBSŁUGI	20
NÁVOD NA OBSLUHU	22
NÁVOD K POUŽITÍ	24
ИНСТРУКЦИ ПО ЭКСПЛУАТАЦИ	26
KULLANMA KILAVUZU	28
NAUDOJIMO INSTRUKCIJA	30
УПЪТВАНЕ ЗА УПОТРЕБА	32
دليل الاستخدام	34

SHAPE SENSE CONTROL 100

013198/A 08/2016 JFS

L I F E I N B A L A N C E

SOEHNLE

Quality & Design by
LEIFHEIT AG
Leifheitstraße 1
56377 Nassau/Germany

www.soehnle.com

- Nicht geeignet für Personen mit elektronischen Implantaten (Herzschrittmacher).
- It is not suitable for persons with electronic implants (pacemaker, etc.).
- Les personnes qui portent des implants médicaux (p. ex. un stimulateur cardiaque etc.) ne doivent pas utiliser cet appareil.
- Non è adatto per persone con impianti elettronici (pacemaker, ecc.).
- Niet geschikt voor personen met elektronische implantaten (pacemaker, enz.).
- No es apropiado para personas con implantes electrónicos (marcapasos, etc.).
- Não é adequado para pessoas com implantes eletrônicos (pacemakers, etc.).
- Ej lämplig för personer med elektroniska implantat (Pace-maker etc).
- Ikke egnet til folk med elektroniske implantater (pacemakere ect).
- Nie nadaje się dla osób, którym wszczepiono elektroniczne implanty (np. rozruszniki serca).
- Prístroj nie je vhodný pre osoby s elektronickými implantátmi (kardiostimulátorom).
- Není vhodný pro osoby s elektronickými implantáty (kardiostimulátory, atd.).
- Никогда не используйте данный прибор вместе с электронными медицинскими приборами, например, кардиостимуляторами.
- Elektronik parçalar (kalp pili, vs.) taşıyan insanlar için uygun değildir.
- Netinka elektroninius implantus (širdies stimulatorius) turintiems žmonėms.
- Уредът не е подходящ за хора с електронни импланти (сърдечни пейсмейкъри).
- غير مناسب للأشخاص الذين يستخدمون أجهزة إلكترونية مزروعة (منظم ضربات القلب).

Alter/Age	♀					♂				
	-	😊	+	++	😊	-	😊	+	++	😊
16-18	< 15	15-25	25-33	> 33	> 52,5	< 8	8-18	18-24	> 24	> 58,5
18-30	< 20	20-29	29-36	> 36	> 49,5	< 8	8-18	18-24	> 24	> 57
30-40	< 22	22-31	31-38	> 38	> 47,5	< 11	11-20	20-26	> 26	> 56
40-50	< 24	24-33	33-40	> 40	> 46,5	< 13	13-22	22-28	> 28	> 55
50-60	< 26	26-35	35-42	> 42	> 46	< 15	15-24	24-30	> 30	> 54
60+	< 28	28-37	37-47	> 47	> 45	< 17	17-26	26-34	> 34	> 53

• Körperfettanteil, body fat content, taux de graisse, percentuale di grasso, lichaamsvetaandeel, procentaje de la grasa corporal, procentagem de gordura corporal, fedtprocent, andelen kroppsfett, test zsirtartalmát, zawartość tłuszczu w ciele, podíl tělesného tuku, долю внутреннего жира, vücuttaki yağ oranını, Küno riebalų procentinė dalis, Процент на телесната мазнина, نسبة الدهون في الجسم

• Körperwasseranteil, body water content, taux d'eau dans le corps, percentuale acqua corporei, lichaamswateraandeel, procentaje del agua corporal, procentagem de água corporal, vandprocent, andelen kroppsvatten, test vızirtalmát, zawartość wody w ciele, podíl vody v těle, тканевой жидкости, ücuttaki su oranını, K no vandens procentinis, Процент на телесната вода, نسبة الماء في الجسم

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A. Operating elements

1. Minus (▼)
2. Confirm (○)
3. Plus (▲)

Sensor Touch: Touch keys only slightly!

B. Preparation

1. Remove battery strip.

Note: Immediately (as long as „0.0“ is still being displayed) place the body analysis scale on a level surface and wait until the scale powers down automatically.

Then begin data input. Otherwise, it might be possible that the weight of the scale would erroneously be added to the initial weighting process.

The same might occur during regular operation if the scale is held in the hand prior to the weighting process.

2. For all measurements, place scale level on a hard surface (not carpet).
3. Cleaning and Care: clean only using slightly damp cloth. Do not use solvents or abrasives. Never submerge scale in water.
4. A country-specific option to switch from kg/cm to lb/in or st/in is possible by pushing the key on the bottom of the scale to switch units.

Attention! Slipping hazard in case of wet surface.

C. Data input

1. Push ○ key.
2. Use ▲/▼ keys to select memory slot and confirm (○).
3. Then use ▲/▼ keys enter body height. Use ○ to confirm.
4. In the next step enter the age and confirm.
5. Finally select sex and confirm.
6. Now immediately place the body analysis scale on the floor, wait until the screen displays „0.0“ and then step on the scale with bare feet.

The initial measurement is indispensable and serves to store the data for later automatic recognition of persons. If no person steps on the scale, data input must be repeated.

The first body analysis is performed after a person steps on the scale.

This scale is equipped with a Time-Out mode. It will automatically switch off after approx. 40 seconds if no keys are pressed.

Data must be entered again after battery change.

D. Body analysis

1. Place scale on hard, level surface and wait until the scale switches off automatically.
2. Step barefoot onto scale. Please stand still.
3. The body weight will be displayed.

The scale determines the body values (body fat, body water).

Note: Body values are calculated based on the determined weight and the personal data.

The scale shuts off automatically.

E. Notifications

1. Replace batteries
2. Overload: max. 180 kg
3. The measured value is higher (H) or lower (L) than the limit.

	♀	♂
	4 - 60 %	10 - 60 %
	5 - 90 %	5 - 90 %

F. Technical data

Max. 180 kg/100 g

Max. 396 lb/0.2 lb

Max. 28 st/0.2 lb

Battery requirement: 1 x 3V CR 2430

Guarantee

Leifheit AG grants you 3 years guarantee as of the date of purchase (or for order upon receipt of the goods) for the quality product at hand. The guarantee claims must be asserted immediately upon occurrence of the defect within the guarantee period. The guarantee is limited to the quality of the products.

The following are excluded from the guarantee:

- (1) wear caused by proper use or natural causes;
- (2) damage caused by improper use or handling (e.g. impact, bump, fall);
- (3) damage caused by non-compliance with stated operating instructions;
- (4) the rechargeable battery

In case of a guarantee claim, Leifheit will at its own discretion provide either repair of defective parts or replace the product. If a repair is not possible and if an identical product is no longer available within the product range for replacement, you will receive a replacement product that is, insofar as possible, of equal value. In case of a guarantee claim, the purchase price cannot be reimbursed. This guarantee also does not grant compensation rights. This guarantee furthermore does not grant any claims for damages.

To assert the guarantee, please contact the seller where you purchased the product and present the defective product, the purchase receipt (copy) and the guarantee slip.

The guarantee is valid globally.

Your statutory rights, in particular warranty rights, continue to be valid and are not limited by this warranty.

EC Conformity

This device complies with the applicable EC Directive 2014/30/EC. You will find the EU compliance declaration on www.soehnle.com.

Battery disposal

EC Directive 2008/12/EC

Batteries are not a part of your regular household waste. You must return batteries to your municipality's public collection or wherever batteries of the respective type are being sold.

Batteries containing hazardous materials are marked with these symbols:

Pb = Battery contains lead

Cd = Battery contains cadmium

Hg = Battery contains mercury

Disposal of electric and electronic devices EC Directive 2002/96/EC

This product is not to be treated as regular household waste but must be returned to a collection point for recycling electric and electronic devices. Further information is available from your municipality, your municipality's waste disposal services, or the retailer where you purchased your product.

Consumer service

If you have questions or comments, the following points of contact are available:

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